




# AseptoActive



Powerful alkaline, chlorine-free cleaner for milking equipment and milk tanks.

- ▲ Water saving due to no pre-rinse
- ▲ Low consumption
- ▲ Free from chlorine

## APPLICATION

After milking, remove as much of the milk as possible (without a water pre-rinse). Clean with **AseptoActive** at 0.25%, minimum 10 min. at minimum 40°C and maximum 75°C. Rinse with water of drinking quality. MUST be used in daily rotation with **HorolithActive**.



## FEATURES AND BENEFITS

**AseptoActive** is an alkaline chlorine-free product with excellent cleaning properties. It has a low usage concentration of 0.25%.

The **Active** products can be used without pre-rinse, enabling to save water up to 35%. They also prevent the release of chlorine into the environment.

## SPECIFICATIONS

AseptoActive	
Physical form	Yellow liquid
Active substance	Sodium and Potassium hydroxide
pH at 20°C:	13.3 - 13.7
Density at 20°C (g/cm <sup>3</sup> ):	1.47 - 1.51
Storage temperature:	-10°C to +40°C

## LOGISTICS

Packaging	25 kg
Product code	2329730
EAN code	4028162329735

A Material Safety Data Sheet is available on request for professional users.  
A Product Data Sheet is available on request.

**Ecolab Europe GmbH**  
Richtistrasse 7  
CH-8304 Wallisellen  
www.ecolab.com