



AseptoSupra

Chlorinated alkaline disinfectant and cleaner for milking equipment and milk tanks.

- ▲ Strong disinfection
- ▲ Very good cleaning performance

APPLICATION

Pre-clean with water. Clean and disinfect with **AseptoSupra** at 0.5 - 0.7 % (minimum 10 min. at 40 - 50°C). Rinse with water of drinking quality. To be used in daily rotation with **HorolithSupra**.



FEATURES AND BENEFITS

AseptoSupra is an alkaline chlorinated product with very good cleaning and disinfection properties. It is suitable for all kinds of standard milking equipment, including cooling tanks. AseptoSupra is in compliance with the Biocide Product Regulation (No. 528/2012).

The **Supra** products are the widely applicable program providing great hygiene results.

SPECIFICATIONS

AseptoSupra	
Physical Form	Light yellow liquid
Active Substance	Sodium hypochlorite
pH at 20°C:	13.5 - 14.0
Density at 20°C (g/cm ³):	1.2 - 1.24
Storage temperature:	-10°C to +30°C

LOGISTICS

Packaging	13 kg	24 kg	72 kg	250 kg
Product code	2331580	2331590	2334660	2331610
EAN code	4028162331585	4028162331592	4028162334661	4028162331615

A Material Safety Data Sheet is available on request for professional users.
A Product Data Sheet is available on request.

Ecolab Europe GmbH
Richtistrasse 7
CH-8304 Wallisellen
www.ecolab.com