

Derma-Perfect®

Laundry washing disinfecting powder
A must for disinfecting cloths used
for teat cleaning

- ▲ Highly effective even at 40°C
- ▲ Designed to be used in washing machine
- ▲ Effective in both hard and soft water areas

APPLICATION



Depending on the level of dirt and the hardness of the water, use between 20 g and 35 g of Derma-Perfect® per kg of cloths. Run a normal washing cycle between 40°C and 60°C following the level of disinfection required.

A pre-wash isn't mandatory.

Cloths can be dried or left damp for the next use.



PREMIUM FEATURES

Rich in active oxygen and cleaning agents, Derma-Perfect® assures a perfect cleaning of the cloths whatever their initial soiling.

Derma-Perfect® is suitable for any kind of washing machine and keep its tremendous properties whatever the hardness of water.



Derma-Perfect® contains peracetic acid, a notified active substance for the disinfection of textiles.

Used at only 60°C, Derma-Perfect® is effective against all kind of bacteria, yeast, spores of *Clostridium difficile* and viruses.



BENEFITS

Derma-Perfect® is a must for a perfect cleaning and disinfection of any kind of cloths.

Its amazing capacity to be used at 40°C or 60°C allows savings in term of electricity and premature wear of the washing machine and of the cloths in comparison to the usual washing protocol at 90°C.

Free from phosphate, Derma-perfect is environmental friendly.

LOGISTICS

Packaging:	20 kg
Product code:	2326250
EAN code:	4028162326253
Palletization:	60 cans

SPECIFICATIONS

Product name :	Derma-Perfect®
Physical Form	Powder
Active Substance	Peracetic Acid
pH Value	10.1 ± 0.5
Storage Temp.	0°C – 40°C
Storage Cond.	Keep away from frost