

# VelouCid® (VelouCid® D)

Double-performing emulsion product to be applied post-milking by dipping.



- ▲ Powerful disinfection by iodine
- ▲ Advanced skin conditioning
- ▲ Optimum protection, in particular during cold weather



## APPLICATION

**VelouCid®** is a ready-to-use product. After each milking apply **VelouCid®** on the whole teat by use of a dipping cup like Ecolab's EasyDip Cup and do not wipe it. Keep the animals standing or walking for 5 min to ensure biocidal effect.

## FEATURES AND BENEFITS

**VelouCid®** provides a fast disinfection against aerobic bacteria (gram positive and gram negative), yeasts, and enveloped viruses according to PTO3 - Veterinary Hygiene (Disinfectants).

**VelouCid®** forms a film which is rich in skin conditioning agents due to the EmulsioCare technology. It ensures optimum protection, in particular during cold weather.

**EmulsioCare**, a true oil in water emulsion, enriches and moisturizes the superficial layers of the skin. **EmulsioCare** gives to **VelouCid®** remarkable dermatological qualities and soothing properties during application.



## SPECIFICATIONS

VelouCid®	
Physical form	Liquid
Active substance	Iodine (PVP)
pH at 20°C	3.8 - 5.0
Density at 20°C	0.98 - 1.02
Storage temperature	5°C to +25°C

## LOGISTICS

Packaging	10 kg	20 kg	58 kg	200 kg	1000 kg
Product code	2233470	2276890	2276920	2276940	2276960
EAN code	4028162233476	4028162276893	4028162276923	4028162276947	4028162276961

A Material Safety Data Sheet is available on request for professional users.  
A Product Data Sheet is available on request.

**Ecolab Europe GmbH**  
Richtstrasse 7  
CH-8304 Wallisellen  
www.ecolab.com