



SIMPLE. CLEAN. EFFICIENT.



Des Spray

SURFACE CLEANER DISINFECTANT

POWERED BY
ECOLAB



DESCRIPTION:

A highly effective ready-to-use cleaner disinfectant with a broad disinfection efficacy.

- Provides excellent detergent action
- Rapid broad spectrum antimicrobial performance for improved safety in working environments including kitchens and public spaces
- Maximum cleaning and hygiene performance



USE INSTRUCTIONS:

1. Always wear PPE when handling the product
2. For application in clean conditions, if visible soiling/debris is present pre-cleaning is a prerequisite, please follow hard surface cleaning procedure
3. Apply ready-to-use solution to surface
4. Wipe surface
5. Leave in contact for the required time shown
6. Rinse food contact surfaces with fresh water and a clean cloth and allow surface to air dry



EFFICACY:

- Bactericidal acc. to EN 1276 dirty conditions for 1 minute
- Bactericidal acc. to EN 13697 dirty conditions for 5 minutes
- Bactericidal & Yeasticidal acc. to EN 16615 dirty conditions for 1 minute
- Yeasticidal acc. to EN 1650 dirty conditions for 1 minute
- Yeasticidal acc. to EN 13697 clean conditions for 1 minute
- Virucidal acc. to EN 14476 enveloped viruses clean conditions for 1 minute
- Virucidal acc. to EN 14476 enveloped viruses dirty conditions for 5 minutes



APPLICATION:

Our certified disinfectant products are 99.99% effective against bacteria and enveloped viruses found in foodservice operations and public spaces when used with the correct cleaning procedure.

For water resistant hard surfaces



SAFETY & STORAGE:

- Do not mix with any other product. Use cleaner disinfectant safely. Always read the label and product information before use. Safety data sheet available for professional users on request.
- Dispose of packaging only if emptied and closed. Product residues: see Safety Data Sheet.
- Active ingredients: 100 g liquid contains: 0.24g Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-C16)) and 0.03g Didecylidimethylammonium chloride (DDAC).
- Store between 0 °C and 40 °C in original closed container only.

PRODUCT CODE: PACK SIZE:

9104560 6x750ml



ADDITIONAL RESOURCES



- MSDS
- WALLCHARTS
- PRODUCT GUIDES



Worldwide Headquarter
Ecolab Inc.
1 Ecolab Place, St. Paul, MN55102
1.800.35.CLEAN

European Headquarter
Ecolab Europe GmbH
Richtistr. 7
CH-8304 Wallisellen